


- Zahnzangen, amerikanische form
- Extracting forceps, american pattern
- Forceps de extracción, forma americana
- Daviers, forme americaine
- Pinze Per Estrazione, forma americana




PI-600-0W1
Fig. W1 winter

obere Schneide- und Eckzähne
upper incisors and canines
incisivos y caninos superiores
incisives et canines supérieures
incisivi e canini superiori




PI-600-001
Fig. 1

obere Schneidezähne und Wurzeln
upper incisors and roots
incisivos y raíces superiores
incisives et racines supérieures
incisivi e radici superiori




PI-600-002
Fig. 2

obere und untere Wurzeln
upper and lower roots
raíces superiores y inferiores
racines supérieures et inférieures
radici superiori e inferiori





PI-600-006
Fig. 6

untere Molaren, Spaltzange
lower molars, splitting forceps
molares inferiores, forceps para separar raíces
molares inferiores, daviers pour separer
des racines
molari superiori, tenaglie per scindere le radici



PI-600-004
Fig. 4

untere und obere Prämolaren
upper and lower premolars
premolares superiores y inferiores
premolaires supérieures et inférieures
premolari superiori e inferiori



PI-600-005
Fig. 5


untere Weisheitszähne
lower third molars
cordales inferiores
denti de sagesse inférieures
denti del giudizio inferiori

- Zahnzangen, amerikanische form
- Extracting forceps, american pattern
- Forceps de extracción, forma americana
- Daviers, forme americaine
- Pinze Per Estrazione, forma americana




PI-600-007
Fig. 7

obere Wurzeln
upper roots
raíces superiores
racines supérieures
radici superiori



PI-600-008
Fig. 8

obere und untere Molaren
upper and lower molars
molares superiores y inferiores
molaires supérieures et inférieures
molari superiori e inferiori




PI-600-013
Fig. 13

obere Schneidezähne
upper incisors
incisivos superiores
incisives supérieures
incisivi superiori



PI-600-010H
Fig. 10H

obere Molaren
upper molars
molares superiores
molaires supérieures
molari superiori



PI-600-010S
Fig. 10S

obere Molaren
upper molars
molares superiores
molaires supérieures
molari superiori



PI-600-014
Fig. 14

untere Schneidezähne
lower incisors
incisivos inferiores
incisives inférieures
incisivi inferiori

- Zahnzangen, amerikanische form
- Extracting forceps, american pattern
- Forceps de extracción, forma americana
- Daviers, forme americaine
- Pinze Per Estrazione, forma americana




PI-600-016
Fig. 16

untere Molaren
lower molars
molares inferiores
molaires inférieures
molari inferiori



PI-600-016S
Fig. 16S

untere Molaren, für Kinder
lower molars , for children
molares inferiores, para niños
molaires inférieures, pour enfants
molari inferiori, per bambini



PI-600-017
Fig. 17

untere Molaren
lower molars
molares inferiores
molaires inférieures
molari inferiori



PI-600-018R
Fig. 18R


obere Molaren, rechts
upper molars, right
molares superiores, derechos
molaires supérieures, droites
molari superiori, destri



PI-600-018L
Fig. 18L

obere Molaren, links
upper molars, left
molares superiores, izquierdos
molaires supérieures, gauches
molari superiori, sinistri

- Zahnzangen, amerikanische form
- Extracting forceps, american pattern
- Forceps de extracción, forma americana
- Daviers, forme americaine
- Pinze Per Estrazione, forma americana



PI-600-069
Fig. 69
Tomes

obere und untere Wurzein
upper and lower roots
raíces superiores y inferiores
racines supérieures et inférieures
radici superiori e inferiori

PI-600-085
Fig. 85

untere Eckzähne, Prämolaren, Molaren
lower canines premolars and molars
caninos, premolares y molares inferiores
canines, prémolaires et molaires inférieures
canini, premolari e molari inferiori

PI-600-085A
Fig. 85A

untere Eckzähne, Prämolaren, Molaren
lower canines premolars and molars
caninos, premolares y molares inferiores
canines, prémolaires et molaires inférieures
canini, premolari e molari inferiori

PI-600-088R
Fig. 88R

obere Molaren, rechts
upper molars, right
molares superiores, derechos
molaires supérieures, droites
molari superiori, destri

PI-600-088L
Fig. 88L

obere Molaren, links
upper molars, left
molares superiores, izquierdos
molaires supérieures, gauches
molari superiori, sinistri


- Zahnzangen, amerikanische form
- Extracting forceps, american pattern
- Forceps de extracción, forma americana
- Daviers, forme americaine
- Pinze Per Estrazione, forma americana



PI-600-150
Cohen
universal
universal
universal
universel
universale




PI-600-099A
Fig. 99A
obere Schneide- und Eckzähne, Prämolaren
upper incisors and canines, premolars
incisivos y caninos, premolares superiores
incisives et canines, prémolaires supérieures
incisivi e canini, premolari superiori



PI-600-099C
Fig. 99C

obere Schneide- und Eckzähne, Prämolaren
upper incisors and canines, premolars
incisivos y caninos, premolares superiores
incisives et canines, prémolaires supérieures
incisivi e canini, premolari superiori



PI-600-151A
Fig. 151A
Cryer

untere Schneidezähne, Prämolaren, Wurzein
lower incisors, premolars, roots
incisivos, premolares, raíces inferiores
incisives, prémolaires, racines inférieures
incisivi, premolari, radici inferiori



PI-600-151
Fig. 151
Cryer

untere Schneidezähne, Prämolaren, Wurzein
lower incisors, premolars, roots
incisivos, premolares, raíces inferiores
incisives, prémolaires, racines inférieures
incisivi, premolari, radici inferiori