

# PASSAU


## IMPEX




# ORTHOPEDIC IMPLANTS



PIO-696



PIO-698  
Basic Fragment Set





## ***ABOUT PRODUCTS***

We are manufacturers of Orthopedic Implants and have vast experience in this line and would like to serve the human beings for this noble cause.

The Orthopedic Implants are produce in a high grade Stainless Steel 316L Titanium Ti6AL4V. The steel have been approved by Industrial Research & Development Centre and also The Pakistan Council of Scientific & Industrial Research Institutions.

These types of Stainless Steel meet all requirements and approved by the ISO committee ISO/TC 150 in the Specification ISO 5832- Part 1,

The Stainless Steel is used for Orthopedic Implants are specified by the National Bureau of Standards in England, united State and Switzerland

Orthopedic Implants has surface finish by the electrically polished and chemecally passivation for remove any substances. During the process in production every Implants checked for accuracy of their gradation, sizes finishing.

### **Analysis of Stainless Steel 316L**

**Chemical composition of Low Carbon Stainless Steel used in Manufacturing of Orthopaedic Implants:**

<b>Carbon</b>	<b>0.03 max</b>	<b>Chromiam</b>	<b>17.00 - 20.00</b>
<b>Manganese</b>	<b>3.00 max</b>	<b>Nickle</b>	<b>10.00 - 14.00</b>
<b>Phosphorous</b>	<b>0.03 max</b>	<b>Molybdenum</b>	<b>2.00 - 4.00</b>
<b>Silicon</b>	<b>0.75 max</b>	<b>Iron</b>	<b>Balance</b>





## CORTICAL SCREW 3.5 mm

1.25 mm 1.75 mm PITCH ARE AVAILABLE

*Same Article Available in Titanium*

Article #	Length	Article #	Length
PIO-700	12 mm	PIO-720	32 mm
PIO-702	14 mm	PIO-722	34 mm
PIO-704	16 mm	PIO-724	36 mm
PIO-706	18 mm	PIO-726	38 mm
PIO-708	20 mm	PIO-728	40 mm
PIO-710	22 mm		
PIO-712	24 mm		
PIO-714	26 mm		
PIO-716	28 mm		
PIO-718	30 mm		



## CORTICAL SCREW 4.5 mm

1.25 mm 1.75 mm PITCH ARE AVAILABLE

*Same Article Available in Titanium*

Article #	Length	Article #	Length
PIO-730	12 mm	PIO-754	36 mm
PIO-732	14 mm	PIO-756	38 mm
PIO-734	16 mm	PIO-758	40 mm
PIO-736	18 mm	PIO-760	42 mm
PIO-738	20 mm	PIO-762	44 mm
PIO-740	22 mm	PIO-764	46 mm
PIO-742	24 mm	PIO-766	48 mm
PIO-744	26 mm	PIO-768	50 mm
PIO-746	28 mm	PIO-770	55 mm
PIO-748	30 mm	PIO-772	60 mm
PIO-750	32 mm	PIO-774	65 mm
PIO-752	34 mm	PIO-776	70 mm



## MALLEOLAR SCREW 4.5 mm

*Same Article Available in Titanium*

Article #	Length
PIO-778	25 mm
PIO-780	30 mm
PIO-782	35 mm
PIO-784	40 mm
PIO-786	45 mm
PIO-788	50 mm
PIO-790	55 mm
PIO-792	60 mm
PIO-794	65 mm
PIO-796	70 mm





## CANCELLOUS SCREW 4.0 mm

Same Article Available in Titanium

Article #	Length	Article #	Length
PIO-798	14 mm	PIO-812	28 mm
PIO-800	16 mm	PIO-814	30 mm
PIO-802	18 mm	PIO-816	35 mm
PIO-804	20 mm	PIO-818	40 mm
PIO-806	22 mm	PIO-820	45 mm
PIO-808	24 mm	PIO-822	50 mm
PIO-810	26 mm	PIO-824	60 mm



## CANCELLOUS SCREW 6.5 mm X 16 mm Thread

Same Article Available in Titanium

Article #	Length	Article #	Length
PIO-826	25 mm	PIO-844	70 mm
PIO-828	30 mm	PIO-846	75 mm
PIO-830	35 mm	PIO-848	80 mm
PIO-832	40 mm	PIO-850	85 mm
PIO-834	45 mm	PIO-852	90 mm
PIO-836	50 mm	PIO-854	95 mm
PIO-838	55 mm	PIO-856	100 mm
PIO-840	60 mm	PIO-858	105 mm
PIO-842	65 mm	PIO-860	110 mm



## CANCELLOUS SCREW 6.5 mm X 32 mm Thread

Same Article Available in Titanium

Article #	Length	Article #	Length
PIO-862	45 mm	PIO-876	80 mm
PIO-864	50 mm	PIO-878	85 mm
PIO-866	55 mm	PIO-880	90 mm
PIO-868	60 mm	PIO-882	95 mm
PIO-870	65 mm	PIO-884	100 mm
PIO-872	70 mm	PIO-886	105 mm
PIO-874	75 mm	PIO-888	110 mm







## CANCELLOUS SCREW 6.5 mm Full Thread

Same Article Available in Titanium



Article #	Length	Article #	Length
PIO-890	25 mm	PIO-910	75 mm
PIO-892	30 mm	PIO-912	80 mm
PIO-894	35 mm	PIO-914	85 mm
PIO-896	40 mm	PIO-916	90 mm
PIO-898	45 mm	PIO-918	95 mm
PIO-900	50 mm	PIO-920	100 mm
PIO-902	55 mm	PIO-922	105 mm
PIO-904	60 mm	PIO-924	110 mm
PIO-906	65 mm		
PIO-908	70 mm		

## LAG SCREW

Same Article Available in Titanium



Article #	Length
PIO-926	50 mm
PIO-928	55 mm
PIO-930	60 mm
PIO-932	65 mm
PIO-934	70 mm
PIO-936	75 mm
PIO-938	80 mm
PIO-940	85 mm
PIO-942	90 mm
PIO-944	95 mm
PIO-946	100 mm
PIO-948	105 mm
PIO-950	110 mm

## WASHERS FOR CANCELLOUS AND CORTICAL SCREW

Article # :PIO-952  
Dia 40 mm  
For Cancellous Screw



Article # :PIO-954  
Dia 65 mm  
For Cancellous Screw & Canulated Screw



## CANULATED SCREW 7.5 mm x 32 mm THREAD

Same Article Available in Titanium

Article #	Length
PIO-956	55 mm
PIO-958	60 mm
PIO-960	65 mm
PIO-962	70 mm
PIO-964	75 mm
PIO-966	80 mm
PIO-968	85 mm
PIO-970	90 mm
PIO-972	95 mm
PIO-974	100 mm
PIO-976	105 mm
PIO-978	110 mm

## CANULATED SCREW 7.5 mm x 16 mm THREAD

Same Article Available in Titanium

Article #	Length
PIO-980	55 mm
PIO-982	60 mm
PIO-984	65 mm
PIO-986	70 mm
PIO-988	75 mm
PIO-990	80 mm
PIO-992	85 mm
PIO-994	90 mm
PIO-996	95 mm
PIO-998	100 mm
PIO-1002	105 mm
PIO-1004	110 mm





**DHS PLATE 95°**  
Same Article Available in Titanium

Article #	Holes	Length
PIO-1006	4 hole	70 mm
PIO-1008	5hole	85 mm
PIO-1010	6 hole	100 mm
PIO-1012	7 hole	115 mm
PIO-1014	8 hole	130 mm
PIO-1016	10 hole	163 mm
PIO-1018	12 hole	198 mm
PIO-1020	14 hole	225 mm
PIO-1022	16 hole	252 mm

**JEWETT NAIL PLATE 120°**  
Same Article Available in Titanium

Article #	Holes	Length	Article #	Holes	Length
PIO-1024	4 hole	70 mm	PIO-1054	6 hole	85 mm
PIO-1026	4 hole	75 mm	PIO-1056	6 hole	90 mm
PIO-1028	4 hole	80 mm	PIO-1058	6 hole	95 mm
PIO-1030	4 hole	85 mm	PIO-1060	8 hole	70 mm
PIO-1032	4 hole	90 mm	PIO-1062	8 hole	75 mm
PIO-1034	4 hole	95 mm	PIO-1064	8 hole	80 mm
PIO-1036	5 hole	70 mm	PIO-1066	8 hole	85 mm
PIO-1038	5 hole	75 mm	PIO-1068	8 hole	90 mm
PIO-1040	5 hole	80 mm	PIO-1070	8 hole	95 mm
PIO-1042	5 hole	85 mm	PIO-1072	10 hole	70 mm
PIO-1044	5 hole	90 mm	PIO-1074	10 hole	75 mm
PIO-1046	5 hole	95 mm	PIO-1076	10 hole	80 mm
PIO-1048	6 hole	70 mm	PIO-1078	10 hole	85 mm
PIO-1050	6 hole	75 mm	PIO-1080	10 hole	90 mm
PIO-1052	6 hole	80 mm	PIO-1082	10 hole	95 mm



**AUSTIN MOORE FORMAL HIP PROSTHESIS**

**NARROW STEM**

Article #	Length	Article #	Length
PIO-1084	38 mm	PIO-1104	48 mm
PIO-1086	39 mm	PIO-1106	49 mm
PIO-1088	40 mm	PIO-1108	50 mm
PIO-1090	41 mm	PIO-1110	51 mm
PIO-1092	42 mm	PIO-1112	52 mm
PIO-1094	43 mm	PIO-1114	53 mm
PIO-1096	44 mm	PIO-1116	54 mm
PIO-1098	45 mm	PIO-1118	55 mm
PIO-1100	46 mm	PIO-1120	56 mm
PIO-1102	47 mm	PIO-1122	57 mm

**REGULAR STEM**

Article #	Length	Article #	Length
PIO-1124	38 mm	PIO-1144	48 mm
PIO-1126	39 mm	PIO-1146	49 mm
PIO-1128	40 mm	PIO-1148	50 mm
PIO-1130	41 mm	PIO-1150	51 mm
PIO-1132	42 mm	PIO-1152	52 mm
PIO-1134	43 mm	PIO-1154	53 mm
PIO-1136	44 mm	PIO-1156	54 mm
PIO-1138	45 mm	PIO-1158	55 mm
PIO-1140	46 mm	PIO-1160	56 mm
PIO-1142	47 mm	PIO-1162	57 mm





# DHS PLATES

www.passauimpex.com

## SEMI TUBULAR PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1164	3 hole	55 mm
PIO-1168	4 hole	71 mm
PIO-1170	5 hole	87 mm
PIO-1172	6 hole	103 mm
PIO-1174	7 hole	119 mm
PIO-1176	8 hole	135 mm
PIO-1178	9 hole	151 mm
PIO-1180	10 hole	167 mm
PIO-1182	12 hole	199 mm



## MICRO PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1184	4 hole	23 mm
PIO-1186	5 hole	29 mm
PIO-1188	6 hole	35 mm



Long Barrel



Short Barrel



## DHS PLATE 135°

Same Article Available in Titanium

Article #	Holes	Length
PIO-1190	4 hole	80 mm
PIO-1192	5 hole	96 mm
PIO-1194	6 hole	112 mm
PIO-1196	8 hole	144 mm
PIO-1198	10 hole	176 mm
PIO-1200	12 hole	208 mm
PIO-1202	14 hole	240 mm
PIO-1204	16 hole	272 mm

135°  
145°  
150°

Long & Short Barrel are Available in all Degrees



## DCP BROAD

Same Article Available in Titanium

Article #	Holes	Length
PIO-1206	5 hole	87 mm
PIO-1208	6 hole	103 mm
PIO-1210	7 hole	119 mm
PIO-1212	8 hole	135 mm
PIO-1214	10 hole	167 mm
PIO-1216	12 hole	199 mm
PIO-1218	14 hole	231 mm
PIO-1220	16 hole	263 mm



## DCP NERROW

Same Article Available in Titanium

Article #	Holes	Length
PIO-1222	4 hole	71 mm
PIO-1224	5 hole	83 mm
PIO-1226	6 hole	103 mm
PIO-1228	7 hole	119 mm
PIO-1230	9 hole	135 mm
PIO-1232	10 hole	167 mm
PIO-1234	12 hole	199 mm
PIO-1236	14 hole	231 mm



## DCP SMALL FREGMENT

Same Article Available in Titanium

Article #	Holes	Length
PIO-1238	5 hole	61 mm
PIO-1240	6 hole	73 mm
PIO-1242	7 hole	85 mm
PIO-1244	8 hole	97 mm
PIO-1246	10 hole	121 mm
PIO-1248	12 hole	145 mm
PIO-1250	14 hole	157 mm



## 1/3 TUBULAR PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1252	3 hole	37 mm
PIO-1254	4 hole	49 mm
PIO-1256	5 hole	61 mm
PIO-1258	6 hole	73 mm
PIO-1260	7 hole	85 mm
PIO-1262	8 hole	97 mm
PIO-1264	9 hole	109 mm
PIO-1266	10 hole	121 mm





## CONDYLAR BLADE PLATE 95°

Same Article Available in Titanium

Article #	Holes	Length
PIO-1268	5 hole	50 mm
PIO-1270	5 hole	60 mm
PIO-1272	5 hole	70 mm
PIO-1274	5 hole	80 mm
PIO-1276	7 hole	50 mm
PIO-1278	7 hole	60 mm
PIO-1280	7 hole	70 mm
PIO-1282		
PIO-1284	9 hole	50 mm
PIO-1286	9 hole	60 mm
PIO-1288	9 hole	70 mm
PIO-1290	9 hole	80 mm
PIO-1292	12 hole	50 mm
PIO-1294	12 hole	60 mm
PIO-1296	12 hole	70 mm
PIO-1298	12 hole	80 mm
PIO-1300	14 hole	50 mm
PIO-1302	14 hole	60 mm
PIO-1304	14 hole	70 mm
PIO-1306	14 hole	80 mm

## CONDYLAR BLADE PLATE 130°

Same Article Available in Titanium

Article #	Holes	Length
PIO-1308	5 hole	50 mm
PIO-1310	5 hole	60 mm
PIO-1312	5 hole	70 mm
PIO-1314	5 hole	80 mm
PIO-1316	7 hole	50 mm
PIO-1318	7 hole	60 mm
PIO-1320	7 hole	70 mm
PIO-1322	7 hole	80 mm
PIO-1324	9 hole	50 mm
PIO-1326	9 hole	60 mm
PIO-1328	9 hole	70 mm
PIO-1330	9 hole	80 mm
PIO-1332	12 hole	50 mm
PIO-1334	12 hole	60 mm
PIO-1336	12 hole	70 mm
PIO-1338	12 hole	80 mm
PIO-1340	14 hole	50 mm
PIO-1342	14 hole	60 mm
PIO-1344	14 hole	70 mm
PIO-1346	14 hole	80 mm





## CASPER PLATE & PEDICAL SCREW



### PEDICAL SCREW

Same Article Available in Titanium

Article #	Length
PIO-1348	30 mm
PIO-1350	35 mm
PIO-1352	40 mm
PIO-1354	45 mm




Same Article Available in Titanium

Article #	Length
PIO-1356	28 mm
PIO-1358	30 mm
PIO-1360	32 mm
PIO-1362	37 mm
PIO-1364	40 mm
PIO-1366	43 mm
PIO-1368	46 mm
PIO-1370	48 mm
PIO-1372	50 mm
PIO-1374	54 mm
PIO-1376	60 mm
PIO-1378	66 mm
PIO-1380	72 mm

### SCHANZ SCREW

Same Article Available in Titanium

Article #	Size	Length
PIO-1382	3.0 mm	90 to 150 mm
PIO-1384	3.5 mm	90 to 150 mm
PIO-1386	4.0 mm	90 to 150 mm
PIO-1388	4.5 mm	90 to 150 mm
PIO-1390	5.0 mm	90 to 150 mm



RUSH NAIL

STEMIEN PIN

SCHANZ SCREW

### STEMIEN PIN

Same Article Available in Titanium

Article #	Size	Length
PIO-1392	3.0 mm	90 to 150 mm
PIO-1394	3.5 mm	90 to 150 mm
PIO-1396	4.0 mm	90 to 150 mm
PIO-1398	4.5 mm	90 to 150 mm
PIO-1400	5.0 mm	90 to 150 mm

### RUSH NAIL

Same Article Available in Titanium

Article #	Size	Length
PIO-1402	3.0 mm	90 to 150 mm
PIO-1404	3.5 mm	90 to 150 mm
PIO-1406	4.0 mm	90 to 150 mm
PIO-1408	4.5 mm	90 to 150 mm
PIO-1410	5.0 mm	90 to 150 mm





Same Article Available in Titanium



PIO-1412 2.0 mm



PIO-1414 2.7 mm



PIO-1416 3.5 mm



PIO-1418 4.5 mm

## RECONSTRUCTION PLATES 4.5 mm Straight

Do not bend more than 15° Fixation with 4.5 & 2.7, 3.5, 2 mm are also Available

## LIMITED CONTACT DYNAMIC COMPRESSION PLATE (LC-DCP) FOR 4.5 mm SCREW ( BROAD)

Fixation with 4.5 & 6.5 mm Screws. Undercuts: For reduced vascular damage. The end holes can be used for cancellous bone screws. The LC-DCP holes allow compression in both directions of the longitudinal axis. Indication: As a tension band plate on the femur and for Humeral Pseudarthrosis (not intended for the Tibia)

Thickness: 6 mm, width: 17.5 mm, Hole Spacing: 18 mm



PIO-1420



PIO-1422

## K-NAIL

PIO-1424 to PIO-1478 A



PIO-1424 to PIO-1478 B




PIO-1424 to PIO-1478 C



## INTERLOCKING NAIL Also Available in Titanium

PIO-1424 to PIO-1478 D



PIO-1424 to PIO-1478 E



Article #	Dia	Length	Article #	Dia	Length	Article #	Dia	Length
PIO-1424	9 mm	34 cm	PIO-1448	10 mm	44 cm	PIO-1472	12 mm	40 cm
PIO-1426	9 mm	36 cm	PIO-1450	10 mm	46 cm	PIO-1474	12 mm	42 cm
PIO-1428	9 mm	38 cm	PIO-1452	11 mm	34 cm	PIO-1476	12 mm	44 cm
PIO-1430	9 mm	40 cm	PIO-1454	11 mm	36 cm	PIO-1478	12 mm	46 cm
PIO-1432	9 mm	42 cm	PIO-1456	11 mm	38 cm			
PIO-1434	9 mm	44 cm	PIO-1458	11 mm	40 cm			
PIO-1436	9 mm	46 cm	PIO-1460	11 mm	42 cm			
PIO-1438	10 mm	34 cm	PIO-1462	11 mm	44 cm			
PIO-1440	10 mm	36 cm	PIO-1464	11 mm	46 cm			
PIO-1442	10 mm	38 cm	PIO-1466	12 mm	34 cm			
PIO-1444	10 mm	40 cm	PIO-1468	12 mm	36 cm			
PIO-1446	10 mm	42 cm	PIO-1470	12 mm	38 cm			



## SMALL T PLATE OBLIQUE ANGLE 3 HOLE HEAD

Same Article Available in Titanium


Article #	Holes	Length
PIO-1480	3 hole	52 mm
PIO-1482	4 hole	63mm
PIO-1484	5 hole	74mm



## SMALL BUTTRESS PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1486	4 hole	54 mm
PIO-1488	5 hole	66 mm
PIO-1490	6 hole	78 mm
PIO-1492	4 hole	54 mm
PIO-1494	5 hole	66 mm
PIO-1496	6 hole	78 mm



## L BUTTRESS PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1498	4 hole	84 mm
PIO-1500	5 hole	101 mm
PIO-1502	6 hole	117 mm
PIO-1504	4 hole	84 mm
PIO-1506	5 hole	101 mm
PIO-1508	6 hole	117 mm

## Y PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1510	3 hole	77 mm
PIO-1512	5 hole	106 mm
PIO-1514	3 hole	77 mm
PIO-1516	5 hole	106 mm







## SMALL T PLATE 3 HOLE HEAD

Same Article Available in Titanium

Article #	Holes	Length
PIO-1518	3 hole	44 mm
PIO-1520	4 hole	56 mm
PIO-1522	5 hole	68 mm
PIO-1524	6 hole	80 mm



## SMALL T PLATE 4 HOLE HEAD

Same Article Available in Titanium

Article #	Holes	Length
PIO-1526	3 hole	44 mm
PIO-1528	4 hole	56 mm
PIO-1530	5 hole	68 mm
PIO-1532	6 hole	80 mm



## T BUTTRESS PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1534	3 hole	81 mm
PIO-1536	4 hole	98 mm
PIO-1538	5 hole	115 mm
PIO-1540	6 hole	132 mm
PIO-1542	7 hole	149 mm
PIO-1544	8 hole	166 mm
PIO-1546	10 hole	200 mm



## T STANDARD PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1548	3 hole	81 mm
PIO-1550	4 hole	98 mm
PIO-1552	5 hole	115 mm
PIO-1554	6 hole	132 mm
PIO-1556	7 hole	149 mm
PIO-1558	8 hole	166 mm
PIO-1560	10 hole	200 mm



## CLOVERLEAF PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1562	3 hole	88 mm
PIO-1564	4 hole	104 mm
PIO-1566	5 hole	120 mm
PIO-1568	6 hole	136 mm



## SPOIN PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1570	5 hole	100 mm
PIO-1572	6 hole	120 mm



## COBRA HEAD PLATE

Same Article Available in Titanium

Article #	Holes	Length
PIO-1574	7 hole	157 mm
PIO-1576	8 hole	170 mm
PIO-1578	9 hole	186 mm
PIO-1580	10 hole	202 mm
PIO-1582	11 hole	218 mm
PIO-1584	12 hole	234 mm

## PLATES FOR DISTAL FEMUR LATERAL

For Use with 4.5 mm Screws

Condylar Buttress Plates

Primary Indication: As a buttress for femoral condyles with multiple fragments ( Left and Right Plates)

Same Article Available in Titanium

Article #	Length	Holes in shaft
PIO-1586	158 mm	7 Right
PIO-1588	158 mm	7 Left
PIO-1590	190 mm	9 Right
PIO-1592	190 mm	9 Left
PIO-1594	221 mm	11 Right
PIO-1596	221 mm	11 Left
PIO-1598	253 mm	13 Right
PIO-1600	253 mm	13 Left
PIO-1602	285 mm	15 Right
PIO-1604	285 mm	15 Left







# BONE INSTRUMENTS

[www.passauimpex.com](http://www.passauimpex.com)



**PIO-1606**  
Double Drill Guide  
4.5 / 3.2



**PIO-1608**  
Double Drill Guide  
6.5 / 3.2



**PIO-1610**  
Universal Drill Guide  
4.5 / 3.2



**PIO-1612**  
Drill & Tap Sleeve  
4.5 mm



**PIO-1614**  
Drill Sleeve  
4.5 mm



**PIO-1616**  
Drill & Tap Sleeve  
6.5 mm



**PIO-1618**  
Holding Sleeve for Screws, for TCP 4.5 / 5.0



**PIO-1620**  
Screwdriver Hexagonal Large,  $\varnothing$  3.5 mm



**PIO-1622**  
Pointed Drill Guide



**PIO-1624**  
T-Handle with Quick Coupling



**PIO-1626**  
Drill Sleeve 5.0, for Drill Bits  $\varnothing$  4.3 mm



**PIO-1628**  
Drill Sleeve Insert



**PIO-1630**  
Countersink, Length 18 cm



**PIO-1632**  
Sharp Hook



**PIO-1634**  
Drill Bit  $\varnothing$  4.5 mm



**PIO-1636**  
Drill Bit  $\varnothing$  3.2 mm



**PIO-1638**  
Tap for Cortex Screw  $\varnothing$  4.5 mm



**PIO-1640**  
Screwdriver Shaft hexagonal  $\varnothing$  3.5 mm



**PIO-1642**  
Depth Gauge for Screw  $\varnothing$  4.5 to 6.5 mm

# BONE INSTRUMENTS

[www.passauimpex.com](http://www.passauimpex.com)



PIO-1644



PIO-1646  
Lexer



PIO-1648  
Lexer



PIO-1650  
Lexer



PIO-1652  
Screw Driver



PIO-1654  
Screw Driver



PIO-1656  
Screw Driver



PIO-1658  
Screw Driver



PIO-1660



PIO-1662  
Mallets



PIO-1664



PIO-1666  
Handle with Quick  
Coupling





# BONE INSTRUMENTS

[www.passauimpex.com](http://www.passauimpex.com)



**PIO-1668**  
Small Wire Cutter  
16 cm Long



**PIO-1670**  
Wire Cutter Large  
Version 22 cm Long



**PIO-1672**  
Flate Nose Parallel Plier  
16 cm Long




**PIO-1674**  
Vice Grip Self Locking  
18 cm Long




**PIO-1676**  
Wire Bending Plier  
15.5 cm Long



**PIO-1678**  
Bone Holding Forcep



**PIO-1682**  
Reduction Forcep  
with Points



**PIO-1684**  
Lane Bone Holding Forcep  
Ranges: 33,45 cm

**PIO-1680**  
Reduction Forcep  
with Points

**PIO-1686**  
Colclough Forceps  
Ranges : 27 cm up

# BONE INSTRUMENTS

[www.passauimpex.com](http://www.passauimpex.com)



PIO-1688



PIO-1690  
Hand Drill with Adjustable  
Handle and 3 Jaw Chuck



PIO-1692  
Tightening  
wrench only

PIO-1694  
Ralk



PIO-1696  
Bunnell





# PASSAU IMPEX

P.O.Box # 1002  
C 21-23 Nishtar Road,  
Small Industrial Estate, Sialkot-51310 Pakistan.  
Tel: +92-52-3555 682 Fax: +92-52-3550 718  
E-mail: [info@passauimpex.com](mailto:info@passauimpex.com)  
Web: [www.passauimpex.com](http://www.passauimpex.com)